

Special Tube Instruction List

Page (1 of 2)

- | | | | | | |
|------------------------------|------------------------------|--------------------------------|-----------------------------|-----------------------------|-----------------------------|
| 1 | 2 | 3 | 4 | 5 | 6 |
| Container Code: GRNLI | Container Code: NAGRN | Container Code: NAGRN10 | Container Code: GOLD | Container Code: YELB | Container Code: YELA |

Volume: 4.0 mL

Volume: 2.0 mL

Volume: 10 mL

Volume: 5.0 mL

Solution B	Solution A
Volume: 6.0 mL	Volume: 8.5 mL

Additive:
Lithium Heparin

Additive:
Sodium Heparin

Additive:
Sodium Heparin

Additive:
**Clot activator & gel
 for serum separation**

Additive: **Sodium
 polyanethol sulfonate.
 Acid citrate dextrose**

Additive:
**Trisodium citrate,
 citric acid, dextrose.**

Tube Color: **Green**



Tube Color: **Green
 with a semi clear cap**

Tube Color: **Green**

Tube Color: **Gold**

Tube Color: **Yellow**


Tube Color: **Yellow**



Special Tube Instruction List

Page (2 of 2)

7	8	9	10	11	12
Container Code: BLU	Container Code: DBLU EDTA	Container Code: DBLU	Container Code: OTH.TB	Container Code: LPST	Container Code: GRAY
Volume: 1.8 mL	Volume: 6 mL	Volume: 6 mL	Volume: 1 mL (each)	Volume: 4 mL	Volume: 4 mL
Additive: Buffered Na Citrate (For low volume)	Additive: K₂ EDTA 10.8 mg Trace Element (Plasma)	Additive: No Additive Trace Element (Serum)	Specific Test Only TB QuantiFERON	Additive: K₂E EDTA Gel separator (Plasma)	Additive: Potassium oxalate/ sodium fluoride
Tube Color: Light Blue with clear cover	Tube Color: Royal Blue	Tube Color: Royal Blue	Tube Colors: Red, Lavendar & Gray with white ring inside	Tube Color: Lavender with yellow ring inside	Tube Color: Gray



Special testing instructions
provided in kit.

