

**Request for Supplies**

**Client:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Address:** \_\_\_\_\_

**Person Requesting:** \_\_\_\_\_ **Phone:** \_\_\_\_\_

**SUPPLY REQUESTS MUST BE SUBMITTED IN WRITING – BY FAX or COURIER DELIVERY**

AQOP provides containers & supplies ONLY for the collection & transport of lab specimens being referred to AQOP for testing. Routine vacutainer tubes (lavender, blue, red, SST, etc) cannot be supplied to facilities who utilize these tubes for specimen collection and testing at their facility.

<b>Item (unit)</b>	<b>QTY ORDERED</b>	<b>QTY SHIPPED</b>	<b>Item (unit)</b>	<b>QTY ORDERED</b>	<b>QTY SHIPPED</b>
<b>SPECIMEN COLLECTION:</b>			<b>AQOP FORMS/GUIDE/LABELS:</b>		
Vacutainer Needle – 21g x 1½" (box/100)		Our staff fills out below:	AQOP Requisitions (ea)		Our staff fills out below:
Vacutainer - Holders (bag/25)			Fee Schedule (ea)		
Vacutainer - blue Citrate NA 1.8 ml (ea)			ABN Advance Beneficiary Notice (20 per bundle)		
Vacutainer - dark green NA Heparin 4 ml (ea)					
Vacutainer - dark green NA Heparin 6 ml (ea)					
Vacutainer - gray top Na Fluoride 2 ml (ea)					
Vacutainer - green PST Gel LI Heparin 3 ml (100/flat or ea)			<b>SPECIMEN TRANSPORT:</b>		
Vacutainer - purple K2 EDTA 3 ml (100/flat or ea)			Biohazard Specimen Transport Bags (50/bundle)		
Vacutainer - red top 8.5 ml Clot Activator (100/flat or ea)			Serum Vials-Plastic (Amber) (ea)		
Vacutainer - yellow top 3.5 ml Clot Activator (100/flat or ea)			Serum Vials-Plastic (approx 100/bag)		
Vacutainer - royal blue K2 EDTA 6 ml (ea)			Trace Element Free Transport Tube (ea)		
Vacutainer - royal blue NO additive 6 ml (ea)					
Vacutainer - yellow ACD Solution A 8.5 ml (ea)					
Quantiferon – TB collection kit			<b>URINE &amp; STOOL COLLECTION:</b>		
			Urine 24-hr container (ea)		
<b>CULTURE/SWABS/PROBE/KIT COLLECTION:</b>			Urine C&S Transfer Straw Kit (Grey Tubes) (ea)		
Blood Culture - Aerobic (ea)			Urine Container - Sterile Cup (ea)		
Blood Culture - Anaerobic (ea)			Occult Blood: <b>Hemocult</b> Cards (ea)		
Blood Culture – Pediatrics, low volume (ea)			Stool for Culture ( <b>orange cap</b> ) (ea)		
GC/Chlamydia vag/cervical Swab collect kit (ea)			Stool for Crypto/Giardia ( <b>prototfix-white cap</b> ) (ea)		
GC/Chlamydia Urine collect kit (ea)			Stool for Ova & Parasite ( <b>prototfix-white cap</b> ) (ea)		
E-Swab - Regular flocced swab <b>White Cap</b> (ea)					
E-Swab - Mini-tipped flocced swab <b>Green Cap</b> (ea)			H.pylori breath kit		
Double Swab for rapid strep– <b>Red Cap</b> Liquid Stuarts (ea)			Saliva Cortisol collect kit		
Flocced <b>NP</b> swab/transport media (UTM) (ea) – <b>For Covid/Influenza/RSV by PCR and Respiratory PCR Panel</b>					
<b>Viral Culture</b> Swab/transport media- (ea)			<b>MISCELLANEOUS SUPPLIES</b> : non-listed supplies may only be furnished on a charge basis, Contact Client Service for questions related to if supply will be charged		
<b>HSV by PCR</b> Swab/transport media (UTM) -					
Vaginitis DNA Probe (VDP) Affirm VPIII (collect kit) (ea)					
Fetal Fibronectin (Spec Collection) Kits (ea)					
<b>AQOP Documentation only</b>					
<b>GLUCOSE TOLERANCE BEVERAGE :</b>			Order Filled by: _____		
Glucola 100g--Orange / Fruit punch / Lemon-lime (ea)-circle 1			Date Filled: ____/____/____		
Glucola 50g-- Orange / Fruit punch / Lemon-lime (ea)-circle 1			Sent By: In-town courier Med Express Same Day Other _____		
Glucola 75g-- Orange / Fruit punch / Lemon-lime (ea)-circle 1					
<b>Revised 11/25/23</b>					